




CelMate® CO<sub>2</sub> Incubators

**CelMate®**

**CO<sub>2</sub> Incubators**  
*Cradle for Beautiful Cells*





## Welcome to Esco

*Esco's Vision is to provide enabling technologies for scientific discoveries to make human lives healthier and safer.*

The Esco Lifesciences Group is committed to deliver innovative solutions for the clinical, life sciences, research, industrial, laboratory, pharmaceutical, and IVF community. With the most extensive product line in the industry, Esco have passed a number of international standards and certifications. Esco represents innovation and forward-thinking designs, that are of the highest standard quality since 1978.

**Availability and Accessibility.** Esco has headquarters in Singapore, Indonesia, and Philippines, with manufacturing facilities are located in Asia and Europe. Research and Development (R&D) is conducted worldwide spanning the US, Europe and Asia. Sales, services and marketing subsidiaries are located in 42 major markets including US, UK, Japan, China and India. Esco regional distribution centers are located in Singapore, Malaysia, Thailand, Vietnam, Myanmar, Indonesia, Philippines, Bangladesh, Hong Kong, Taiwan, South Korea, China, Japan, India, UAE, Central and South Africa, Denmark, Germany, Italy, Lithuania, Russia, United Kingdom, and USA. Because of our worldwide presence, you can be sure that Esco is within your reach.

**High Quality, Reliable, and Dependable.** Esco products are of high quality, reliable, and dependable; assuring customers of research accuracy. Cross functional teams from Esco Production, R&D, Quality Assurance, and Senior Management, are regularly assembled to review and implement areas for improvement.

**Esco Cares for Your Safety.** Esco focuses on providing safety not just for your samples but also for you and the environment.

**Esco Cares for Your Comfort.** Building ergonomic designs and reducing noise levels of the units ensures comfort for our users.

**Esco Cares for the Environment.** One in every four of Esco's employees is involved in R&D and a number of them evaluate new components and/or designs to produce energy efficient equipment. Being GREEN is more than just modifying parts used to produce a new energy efficient technology, it is also embodied in the every aspect of the company.

**Customer Service and Support.** Our service does not stop once purchase has been done. Esco gives on-time customer service and offers end-user seminars, service training, preventive maintenance, and provides educational materials and informative videos.

As Esco takes the opportunity to respond to the world's needs, we aim not only to contribute in the advancement of scientific discoveries but also in making the world a safer, healthier, and better place to live in.

## Products and Application

### Laboratory Equipment

#### Sample Handling and Preparation

- Class I Biological Safety Cabinets
- Class II Biological Safety Cabinets
- Class II Type A2 Biological Safety Cabinets
- Class II Type B1 Biological Safety Cabinets
- Class II Type B2 Biological Safety Cabinets
- Class III Biological Safety Cabinets
- Horizontal Laminar Flow Cabinets
- Vertical Laminar Flow Cabinets
- Laboratory Animal Research Workstations
- Laboratory Centrifuges

#### Sample Cultivation

- CO<sub>2</sub> Incubators, Direct Heat Air-Jacketed
- CO<sub>2</sub> Incubators with Cooling System
- CO<sub>2</sub> Incubators with High Heat Sterilization
- Laboratory Shakers

#### Amplification and Detection

- Conventional Thermal Cyclers
- Microplate Shakers
- PCR Cabinets

#### Sample Storage & Sample Protection Solutions

- Ultra-low Temperature Freezers
- Lab Refrigerators and Freezers
- Sample Database Management Software
- Intelligent Remote Monitoring Application Protocol
- Remote Monitoring, Datalogging, Programming Software
- Wireless Monitoring System

#### Chemical Research

- Ducted Fume Hoods
- Ductless Fume Hoods
- Filtered Storage Cabinets
- Powder Weighing Balance Enclosure
- Exhaust Blowers
- Fume Hood Airflow Monitor

#### General Equipment

##### Laboratory Thermostatic Products

- Forced Convection Laboratory Oven
- Forced Convection Laboratory Incubator
- Natural Convection Laboratory Incubator
- Refrigerated Laboratory Incubator

### Medical / IVF Equipment

#### Controlled Embryo Handling

- Esco Multi-Zone ART Workstation
- Esco Multi-Zone ART Workstation Class II
- AVT Anti-Vibration Table
- Semi-Closed Environment (SCE) IVF

#### Safe Embryo Culture

- MIRI® Multiroom Incubator
- MIRI® II Multiroom Incubator
- Mini MIRI® Humidified Incubator
- Mini MIRI® Dry Incubator
- CelCulture® CO<sub>2</sub> Incubator

#### Innovative Time-Lapse Imaging

- MIRI® Time-Lapse Incubator

#### Accurate Quality Control

- MIRI® GA Gas and Temperature Validation Unit

#### Unique Consumables

- CultureCoin®

### Healthcare

#### Esco Pharma Products

##### Airflow Containment

- BioBooth™
- Ceiling Laminar Airflow (CLAF)
- Cytoculture® Cytotoxic Safety Cabinet
- Pharmacon™ Downflow Booth
- Esco Garment Storage Cabinet
- Esco Glassware Hoods
- Laminar Flow Horizontal/Vertical Trolley (LFH/VT)
- Laminar Flow Straddle Units
- Evidence Drying Cabinet

##### Isolation Containment

- Advanced Processing Platform Isolator (APPI)
- Aseptic Containment Isoaltor (ACTI)
- Blood Cell Labelling Isolator
- Streamline® Closed Restricted Access Barrier System (SLC-RABS)
- Containment Barrier Isolator (CBI)
  - CBI-Unidirectional (CBI-U)
  - CBI-Turbulent (CBI-T)
  - CBI-Class III Biosafety Cabinet (CBI-III)
  - CBI-Convertible Class III/Class I Biosafety Cabinet (CBI-H)
- Isoclean® Healthcare Platform Isolator (HPI)
  - HPI-G3-Without Filter Below Work Zone
  - HPI-G3-With Filter Below Work Zone
  - HPI-Inflatable Seal (HPI-IS)
- General Processing Platform Isolator
  - GPPI-Inflatable Seal (GPPI-IS)
  - GPPI-Static Seal (GPPI-SS)
- Streamline® Compounding Isolator
  - SCI - Isolator Configuration
  - SCI - Class III Biosafety Cabinet (SCI-III)
- Technetium Dispensing Isolator
- Turbulent Flow Aseptic Isolator (TFAI)
- Weighing and Dispensing Containment Isolator (WDCI)

##### Cross Contamination Facility Integrated Barrier

- BioPass™ Pass Through
- Cleanroom Air Showers
- Dynamic Pass Boxes/ Dynamic Floor Laminar Hatches
- Infinity® Air Shower Pass Box
- Esco Sputum Booth
- Infinity® Pass Boxes
- Infinity® Cleanroom Transfer Hatch
- Soft Capsule® Soft Wall Cleanroom

##### Ventilation Containment

- Ventilated Balance Enclosure

#### Esco VacxiXcell Products

##### Bioreactors and Fermenters

- CelXrocker™
- CelCradle™
- CelShaker™
- CelCradle™ X
- CelCradle Semi-Automated Harvesting System™ (CCX-SAH)
- BioXcell™
- StirCradle™
- StirCradle™ PRO
- TideXcell™
- TideXcell™ Harvesting System (TXLHS)
- VXL™ Hybrid Bioreactor

##### Cell Culture Monitoring, Media and Consumables

- Super Plus™
- Plus™ Vero
- Plus™ MDCK
- Plus™ MDCK II
- BioNOC™ II macrocarriers
- GlucCell™ Glucose Monitoring System
- CVD Kit

##### Filling Line Equipment

- Filling Line Isolators
- cRabs (close restricted access barriers)
- oRabs (open restricted access barriers)

##### Integrated Solutions

- Cell Processing Isolator
- Cell Processing Center

#### Esco TaPestle Rx Products

##### Pharmacy Compounding Solutions

- Compounding Pharmacy Isolators (SCI, HPI, CBI, GPPI)
- Safety Cabinets and Enclosures (CYT, Class II BSC, VBE, LFC)
- Aseptic Filling Systems

##### Radiopharmacy Equipment Solutions

- Radioisotope Fume Hood
- Lead-lined Biosafety Cabinet
- Technetium Dispensing Isolator
- Blood Cell Labeling Isolator
- GMP-compliant Radioisotope Dispensing Isolator



# CelMate®

## CO<sub>2</sub> Incubators

### INTRODUCTION

Esco now offers the new CelMate®, 50 L, 170 L, and 240 L, entry-level cell culture CO<sub>2</sub> incubator with superb contamination control. It is specifically designed for laboratories looking for a cost-effective CO<sub>2</sub> incubator that can provide the best protection for their cell culture.

Sleek, reliable, and intuitive, Esco CelMate® CO<sub>2</sub> incubators provide all-rounded sample protection that brings your scientific dreams one step closer to reality.


## KEY FEATURES

### CelMate® CO<sub>2</sub> INCUBATORS

*Cradle for Beautiful Cells*



CelMate® CO<sub>2</sub> Incubators available in  
50 L, 170 L, and 240 L



### SHELVING

- Perforated shelving to improve uniformity
- Anti-tip
- Stainless steel
- Built-in grip
- Dismantles without tools for easy cleaning

### DIRECT HEAT & AIR JACKET

- Fast and uniform heating
- Rapid temperature recovery without overshoot
- Air jacket improves chamber stability



### DUCT WORK

- Directs air flow for rapid recovery and excellent uniformity
- Easily removed for cleaning



### WATER PAN

- Precisely heated by base heater to provide high humidity
- Gentle airflow over water surface accelerates humidity recovery



### ROUNDED CORNERS

- Seamless design
- Facilitates easier cleaning



## ULPA FILTER

- 99.999% efficient, superior to conventional HEPA filters
- Filters air continuously
- Chamber returns to ISO Class 5 cleanliness in 11 minutes upon door closing to prevent contamination



## CO<sub>2</sub> SENSOR

- Heat-resistant IR sensor
- Single-beam, dual-wavelength IR sensor is drift-free
- Auto-zeroing

## DOOR SWITCH

Automatically turns off the blower, heater, UV, and gas supply when the door is opened

## TOP COVER

Provides quick access to electrical panel components

## SMARTSENSE™ MICROCONTROLLER INTERFACE

Intuitive, fully equipped control and monitoring system

## BLOWER\*

Gentle airflow in chamber improves recovery and uniformity

*\*Not available in 50 L Model (including top, back, and bottom plenum)*

## OUTER DOOR

- Reversible
- Heated to prevent condensation

## SAMPLE PORT

Allows direct measurement of chamber atmosphere such as CO<sub>2</sub> and O<sub>2</sub> concentration

## INNER GLASS DOOR

For observing sample cells inside the chamber during operation

## DOOR LATCH

To lock / unlock the glass door

## LEVELING FEET


Easily adjustable

## QUALITY ESCO CONSTRUCTION


- Electrogalvanized steel with white oven-baked epoxy-polyester antimicrobial powder-coated finish.
- External surfaces are powder coated with Esco **ISOCIDE™** to eliminate 99.9% of surface bacteria within 24 hours of exposure.
- Ensures a healthier, safer and cleaner lab environment.

# VIVOCELL™ PRECISE PARAMETER CONTROL

## BEST UNIFORMITY AND CONTROL AMONG THE COMPETITION




Different lines represent different sensor positions inside the chamber. Esco CelMate® has uniformity variance of less than  $\pm 0.5^\circ\text{C}$  which means all the samples are evenly heated.\*



Minimal fluctuation ( $\pm 0.2^\circ\text{C}$ ) ensures temperature stability.\*

## FAST CO<sub>2</sub>, TEMPERATURE AND HUMIDITY RECOVERY WITHOUT OVERTHOOT





Precisely-tuned sensor and software result in fast recovery of CO<sub>2</sub> without overshoot. This ensures uniform CO<sub>2</sub> levels even with frequent incubator door opening.

Recovery of both temperature and humidity is twice as fast as conventional incubators.

- Company A's model: overshoot.
- Company B's model: slow recovery.
- Esco CelMate®: fast recovery, no overshoot.

## DIRECT HEAT AND AIR JACKET





- Direct heating enables rapid temperature recovery while air jacket provides isolation against ambient temperature fluctuations.
- Precise heating in the chamber is achieved by using 8 heaters located in 3 zones. The 3 zones are intelligently controlled by the microcontroller for best temperature uniformity and minimal fluctuation.

*Note: For CLM-050B-\_\_ models, heaters are located in 4 zones (main, bottom, outer door, and back).*

- The main heater provides precise temperature control.
- The bottom heater warms the water pan and provides humidity.
- The outer door heater prevents condensation on glass door and facilitates temperature recovery.

## VENTIFLOW™ FORCED CONVECTION




- No disturbance to cell culture.
- Blower automatically stops when door is opened to minimize mixing of chamber and room air.
- Accelerates recovery of chamber air to ISO Class 5 Cleanliness after door closing to prevent contamination.
- Improves CO<sub>2</sub>, humidity and temperature uniformity.
- Filtered air circulates across water pan to accelerate humidifying process.

*\*Units were factory-tested under controlled environmental conditions per DIN 12880 standard and Esco method. Esco does not guarantee identical results in the field under differing conditions. Original report available upon request. Model used in the test is CLM-170B-8.*

# ROBUST CONTAMINATION CONTROL


## STERISAFE™ ULPA FILTRATION SYSTEM



- Chamber air is continuously filtered by ULPA filters to keep the chamber at ISO Class 5 cleanliness. This ensures that all contaminants from both room air and chamber air are filtered, thus only clean air is recirculated.
- ULPA filters operate at 99.999% efficiency, superior to conventional HEPA filters which are 99.99% efficient.
- Chamber achieves ISO Class 5 cleanliness 11 minutes after door closing.\*

\*Units were factory-tested under controlled environmental conditions per DIN 12880 standard and Esco method. Esco does not guarantee identical results in the field under differing conditions. Original report available upon request. Model used in the test is CLM-170B-8.

## VALIDATED SWIFTCON™ OVERNIGHT DECONTAMINATION CYCLE



Microorganisms	Before Decon	After Decon
<i>Bacillus atrophaeus</i>	1.59 x 10 <sup>6</sup>	0
<i>Aspergillus brasiliensis</i>	1.52 x 10 <sup>4</sup>	0
<i>Pseudomonas aeruginosa</i>	2.38 x 10 <sup>6</sup>	0
<i>Staphylococcus epidermis</i>	2.33 x 10 <sup>6</sup>	0
<i>Escherichia coli</i>	1.57 x 10 <sup>6</sup>	0
<i>Staphylococcus aureus</i>	5.72 x 10 <sup>6</sup>	0
<i>Enterobacter faecalis</i>	2.15 x 10 <sup>6</sup>	0

- The automated SwiftCon™ 90°C moist heat decontamination cycle has been proven effective in deactivating normally resistant fungi, bacterial spores and vegetative cells by the Health Protection Agency (HPA) in UK.
- Full decontamination cycle completes within 20 hours.
- Independently proven to be as effective as high temperature decontamination.
- Lower temperature causes less damage to electronic components, therefore prolongs the life span of the incubator.


## GAS INJECTION LINES ARE FILTERED



- All gas injection lines are filtered via 0.2 micron inlet filters to remove impurities and contaminants before being injected into the chamber.
- Inlet filters are field-replaceable and are located external to the incubator.

# CONTROLLER TYPE

## USER - FRIENDLY SOFTWARE INTERFACE



1. **Start / stop decontamination cycle**
2. **HEAT LED**  
Lights up when heat is applied to chamber
3. **°C is lit** when displaying the temperature
4. **% RH is lit** when displaying the humidity level
5. **Enter menu / go back to previous menu**
6. **Scroll up / increase value**
7. **ALARMS LED**  
Will blink when errors and warnings occur
8. **Mute alarms**
9. **INJECT LED**  
Lights up when gas is injected
10. **% O<sub>2</sub> is lit** when displaying the O<sub>2</sub> concentration  
(not applicable to CelMate®)
11. **% CO<sub>2</sub> is lit** when displaying the CO<sub>2</sub> concentration
12. **Confirm value / enter a menu**
13. **Scroll down / decrease value**

### • Comprehensive, user-configurable alarms:

- Temperature
- CO<sub>2</sub>
- Humidity (if installed)

### • CelAlert™ alarm system reminds user to replace parts.

```
NEW CO2 TANK
CO2 TANK NET WEIGHT
IN KG: 0031
NEW CO2 TANK IS SET
```

In addition to CO<sub>2</sub> tank low alarm, CelAlert™ has CO<sub>2</sub> tank depletion reminder that automatically calculates how much CO<sub>2</sub> gas is left in the tank and alerts user to replace the tank one week before the gas is depleted. This gives the user some buffer time to place orders for new tanks.

```
SET ULPA REMINDER
ULPA REMINDER
TIME <MONTHS> :12
REMINDER IS SET
```

ULPA reminder will alert user to replace ULPA filter.

### • Intelligent data and event logger records all incubator parameters for on-screen recall. A 2 MB built-in flash memory guarantees long-term storage of data.

```
MAIN MENU
CALIBRATION
ADMIN
>DATALOGGING
```

```
DATALOGGING
>SHOW TEMP DATA LOG
SHOW %CO2 DATA LOG
DATA LOG PERIOD
```

```
012016 0724 36.8°C
012016 0719 37.0°C
012016 0714 37.1°C
012016 0709 37.3°C
```

### • Diagnostic interface and online quick help provide comprehensive solutions to frequently encountered problems.

### • New Safety Alarm System

- Temperature Fail-Safe System prevents overshooting of temperature display to +0.4°C of the set point.
- The Auto-Reset Watchdog will automatically reset the system in the unlikely event of system failure, preventing the controller from freezing.
- %CO<sub>2</sub> Failure Mode prevents build-up of %CO<sub>2</sub> over set point in cases of CO<sub>2</sub> sensor defect.

```
SYSTEMS SETPOINTS
TEMP=37.0°C
CO2= 5.0%
BACKUP RUNNING
```

```
SYSTEMS SETPOINTS
TEMP=37.0°C
CO2= 5.0%
CHECK CO2 SYSTEM
```



## Voyager

### Remote Monitoring, Datalogging, Programming Software

Esco Voyager is a PC-based software package developed for the remote monitoring, datalogging, and programming / device configuration of Esco thermostatic products.

A centralized monitoring and control system for the laboratory, Esco Voyager provides extra protection for you and your samples.


Voyager interfaces with individual Esco equipment over RS485 using the EscoBUS communications protocol. Multiple equipment may be interfaced to a single PC.

### Compatible Equipment

- CelCulture® - CO<sub>2</sub> Incubator (CCL)
- CelMate® - CO<sub>2</sub> Incubator (CLM)
- Lexicon® - Ultra-low Temperature Freezer
- Isotherm® - Forced Convection Oven (OFA)
- Isotherm® - Forced Convection Incubator (IFA)
- Isotherm® - Refrigerated Incubator (IFC)
- Isotherm® - Natural Convection Incubator (INA)



# REAR PANEL



## 1 Cooling Fan

Prevents the electrical panel from overheating.



## 5 Alarm Contact

A set of relay contacts located on the rear panel of the unit is provided to monitor temperature, humidity, CO<sub>2</sub> alarms. These can be connected to a remote alarm system.



## 2 Power Supply Inlet

Connects the incubator unit to the power source.



## 6 CO<sub>2</sub> Gas Supply Inlet

Connects the CO<sub>2</sub> gas supply to the incubator. Inlet pressure requirement is 15 psi.




## 3 RS485 Communication Port

Provides serial communication port for PC. It can be daisy-chained from one product to another and can also be connected to a PC




## 7 Gas Inlet Filter

Provided to remove any contaminants from the gas supply.



## 4 Analog Port (Optional)

Allows the incubator to output analog signals representing temperature, CO<sub>2</sub>/O<sub>2</sub> concentration and relative humidity, depending on the options available in the incubator. This allows the incubator to be connected to an in-house data acquisition or alarm system.



## 8 Access Port

Allows cables, hoses or additional sensors to be routed into the work space. A rubber stopper with controlled leak is installed as standard configuration and is part of standard accessories.

# ORDERING INFORMATION

## IR SENSOR MODEL WITH STAINLESS STEEL CHAMBER


MODELS	ITEM CODE	DESCRIPTION
CLM-050B-8	2170269	CelMate® Incubator 50 L, IR Sensor, CO <sub>2</sub> Control, Moist Heat Decon, 230 VAC, 50/60 Hz (Without Decon Pump)
CLM-050B-9	2170345	CelMate® Incubator 50 L, IR Sensor, CO <sub>2</sub> Control, Moist Heat Decon, 115 VAC, 50/60 Hz (Without Decon Pump)
CLM-170B-8	2170106	CelMate® Incubator 170 L, IR Sensor, CO <sub>2</sub> Control, ULPA, Moist Heat Decon, 230 VAC, 50/60 Hz (Without Decon Pump)
CLM-170B-9	2170250	CelMate® Incubator 170 L, IR Sensor, CO <sub>2</sub> Control, ULPA, Moist Heat Decon, 115 VAC, 50/60 Hz (Without Decon Pump)
CLM-240B-8	2170107	CelMate® Incubator 240 L, IR Sensor, CO <sub>2</sub> Control, ULPA, Moist Heat Decon, 230 VAC, 50/60 Hz (Without Decon Pump)
CLM-240B-9	2170251	CelMate® Incubator 240 L, IR Sensor, CO <sub>2</sub> Control, ULPA, Moist Heat Decon, 115 VAC, 50/60 Hz (Without Decon Pump)

## IR SENSOR MODEL WITH STAINLESS STEEL CHAMBER AND UV LAMP


MODELS	ITEM CODE	DESCRIPTION
CLM-170B-8-UV	2170252	CelMate® Incubator, 170L, IR sensor, CO <sub>2</sub> control, Moist Heat Decon, UV, 230VAC, 50/60 Hz (Without Decon Pump)
CLM-240B-8-UV	2170254	CelMate® Incubator, 240L, IR sensor, CO <sub>2</sub> control, Moist Heat Decon, UV, 230VAC, 50/60 Hz (Without Decon Pump)
CLM-170B-9-UV	2170253	CelMate® Incubator, 170L, IR sensor, CO <sub>2</sub> control, Moist Heat Decon, UV, 115VAC, 50/60 Hz (Without Decon Pump)
CLM-240B-9-UV	2170255	CelMate® Incubator, 240L, IR sensor, CO <sub>2</sub> control, Moist Heat Decon, UV, 115VAC, 50/60 Hz (Without Decon Pump)

## ENGINEERING DRAWING


### MODEL 50L



### MODEL 170L



### MODEL 240L



1. Control Panel
2. On / Off Switch
3. Small Fan (50L) / Blower (170L & 240L)
4. Adjustable Shelves
5. IR Sensor

6. Access Port
7. Humidity Pan
8. CO<sub>2</sub> Gas Supply
9. Alarm Contact
10. Analog Output

11. RS485
12. Power Supply Inlet
13. Cooling Fan

# GENERAL SPECIFICATIONS

CELMATE® CO<sub>2</sub> INCUBATORS

GENERAL SPECIFICATIONS		CLM-050B-		CLM-170B-		CLM-240B-	
CELMATE® CO <sub>2</sub> INCUBATORS							
TEMPERATURE							
Temperature Control Method		Direct Heat and Air Jacket using Microcontroller PI					
Ambient Temperature Range		18 to 34°C (64 to 93 °F)					
Temperature Range, °C		Ambient +7 to 60					
Temperature Uniformity, °C*		< ±0.5					
Temperature Accuracy, °C*		<± 0.1					
Temperature Recovery Time** (after 30 seconds door opening, 98% from initial value)		≤5 minutes		≤5 minutes		≤6 minutes	
CO <sub>2</sub>							
CO <sub>2</sub> Control System		Microcontroller PI					
CO <sub>2</sub> Range, % CO <sub>2</sub>		0-19.5					
CO <sub>2</sub> Accuracy, % CO <sub>2</sub> ***		±0.1					
CO <sub>2</sub> Sensor		Infrared (IR) Sensor					
CO <sub>2</sub> Recovery Time**** (after 30 seconds door opening, 98% from initial value)		≤6 minutes		≤5 minutes		≤6 minutes	
HUMIDITY							
Humidification Method		Humidity pan					
Humidity Range, %RH (at 37°C)*****		85-90					
PHYSICAL CONSTRUCTION							
Interior Volume		50 L (1.8 ft³)		170 L (6 ft³)		240 L (8.5 ft³)	
External Dimensions (W x D x H)		500 x 500 x 655 mm (19.7" x 19.7" x 25.8")		660 x 660 x 900 mm (26.0" x 26.0" x 35.4")		750 x 770 x 900 mm (29.5" x 30.3" x 35.4")	
Internal Dimensions (W x D x H)		345 x 375 x 388 mm (13.6" x 14.8" x 15.3")		505 x 535 x 633 mm (19.9" x 21.1" x 24.9")		595 x 640 x 633 mm (23.4" x 25.2" x 24.9")	
Chamber Construction	Main Body	Electrogalvanized steel with Isocide™ antimicrobial coating					
	Interior Material	Stainless steel, type 304					
	Number of Shelves	3		4		4	
	Maximum Number of Shelves	4		7		7	
	Shelves Area (W x D)	300 x 335 mm (11.8" x 13.2")		465 x 470 mm (18.3" x 18.5")		550 x 560 mm (21.7" x 22.0")	
	Maximum Load per Shelf	4 kg/shelf (8.8 lbs/shelf)		11 kg/shelf (24.3 lbs/shelf)		15 kg/shelf (33.1 lbs/shelf)	
Electrical Configuration 110-130 VAC, 50/60 Hz	Nominal Power at 37°C	40.9 W		42.2 W		42.2 W	
	Maximum Power Consumption	675.5 W		1184.3 W		1727.9 W	
	Full Load Amps	5.3 A		9.2 A		13.4 A	
Electrical Configuration 220-240 VAC, 50/60 Hz	Nominal Power at 37°C	40.9 W		42.2 W		42.2 W	
	Maximum Power Consumption	598.8 W		1008.9 W		1270 W	
	Full Load Amps	2.5 A		4.2 A		6.5 A	
Net Weight		52 kg (114.6 lbs)		101 kg (222.67 lbs.)		121 kg (266.76 lbs.)	
Shipping Weight		70 kg (154.3 lbs)		120 kg (264.6 lbs)		155 kg (341.7 lbs)	
Shipping Dimensions (W x D x H)		660 x 650 x 900 mm (26.0" x 25.6" x 35.4")		850 x 720 x 1120 mm (33.5" x 28.3" x 44.1")		850 x 850 x 1120 mm (33.5" x 33.5" x 44.1")	
Shipping Volume		0.39 m³ (13.7 ft³)		0.70 m³ (24.85 ft³)		0.79 m³ (28.03 ft³)	
CONTAMINATION CONTROL							
Contamination Control Methods		1) Main body is electro-galvanized steel with Isocide™ antimicrobial coating; 2) 90°C moist heat OVERNIGHT decontamination cycle (HPA validated); 3) ULPA filter***** 4) 0.2-micron inlet filter for gas inputs; 5) 1-micron air circulation filter 6) UV Lamp (for CLM-UV models only)					

All data recorded is specified for standard models with unloaded chambers and tested under optimum factory setting conditions of 23°C and 60% ambient humidity.

\*Results are achieved when tested at 37 °C as set point. Results may vary if set point changes and calibration is needed.

\*\*For temperature not exceeding 37.2 °C.

\*\*\*Results are achieved when tested at 5% CO<sub>2</sub> as set point. Results may vary if set point changes and calibration is needed.

\*\*\*\*For CO<sub>2</sub> level not exceeding 5.2%.

\*\*\*\*\* Esco does not guarantee condensation-free chamber at humidity level higher than 90%.

\*\*\*\*\*Not available for 50 L models.

# OPTIONS AND ACCESSORIES




## COA-1001 / COA-1001-F Humidity Display

This option allows the incubator to monitor the relative humidity inside the chamber. The probe for the sensor works in freezing conditions (-70°C) and also in temperatures up to 180°C. The sensor is easy to install and has excellent accuracy. The airflow in the chamber does not affect the measurement. The sensor is maintenance-free. It does not need to be removed during 90°C moist heat decontamination cycle.



## COA-1002 / COA-1002-F CO<sub>2</sub> Backup


This option allows two tanks of CO<sub>2</sub> to be connected to the incubator. It will automatically switch from the primary tank to the secondary tank when low gas pressure is detected on the primary tank.



## COA-1005 / COA-1005-F Analog Output

A set of relay contacts is provided at the rear of the incubator that allows the incubator to output analog signals representing the temperature, CO<sub>2</sub> / O<sub>2</sub> content and relative humidity, depending on the options available in your incubator. This allows the chamber to be connected to an in-house data acquisition or alarm system. This option can also be field-installed.

The analog signal outputs can be set to operate in either voltage DC (0-5 Vdc) or current (4-20 mA) mode. The factory default setting is voltage. Switch on the board to toggle between the modes.



## COA-2018-F (50L) / COA-2001-F (170 L) / COA-2019-F (240 L) Roller Base

Roller base is available with casters for mobility of your incubators and to provide protection against floor contamination.



## COA-2020-F (50L) / COA-2002-F (170 L) / COA-2021-F (240 L) Floor Stand 200 mm (8.0") With Adjustable Feet

Floor stands are available with adjustable feet, with a nominal range of 180 mm to 250 mm (7.1" to 9.8") for comfortable access to the incubator and to avoid floor contamination.



## COA-2022-F (50L) / COA-2003-F (170 L) / COA-2023-F (240 L) Floor Stand 700 mm (27.6") With Casters

This support stand raises the incubator to a height of 700 mm (27.6") above the floor for comfortable access. It comes with casters for mobility of your incubators.



## COA-2005-F 2-Stage Gas Regulator for CO<sub>2</sub>

CO<sub>2</sub> and N<sub>2</sub> gas input regulators reduce pressure from the tank to the incubator. It has dual pressure gauges, barbed line connection and shut-off valve. It prevents over-pressurization of the gas supply into the incubator which could cause the tubing to burst.

- CGA 320 connector (U.S. Standard)
- BP-B5341-#8-NT4 connector (British Standard)

*Note: Compatible with European DIN477, French NFE29-650 and Australia AS2473*

- G5/8-RH connector (China Standard)





### COA-2024-F (50L)/ COA-2007-F (170 L)/ COA-2025-F (240 L) Extra Shelf (Stainless Steel) for Standard Stainless Steel Chamber

Each CelMate® CO<sub>2</sub> incubator comes standard with 3 shelves for 50 L / 4 shelves for 170 L & 240L and it can accommodate up to a maximum of 4 shelves for 50 L / 7 shelves for 170 L & 240 L.



### COA-2010-F Electronic CO<sub>2</sub> Analyzer, For CO<sub>2</sub> / Temp Measurement

### COA-2016-F Electronic CO<sub>2</sub> + O<sub>2</sub> Analyzer, For CO<sub>2</sub> / O<sub>2</sub> / Temp Measurement

### COA-2017-F Electronic CO<sub>2</sub> + O<sub>2</sub> + RH Analyzer, For CO<sub>2</sub> / O<sub>2</sub> / RH / Temp Measurement

The electronic analyzer allows the measurement of CO<sub>2</sub> concentration, O<sub>2</sub> concentration, relative humidity and temperature (temperature probe already included).



### COA-2012-F 6" Chart Recorder, Temp, 115/230 VAC, 50/60 Hz

The chart recorder provides an easy-to-read graph of data vs time. It is a reliable, accurate, and stable instrument for on-the-spot written documentation of incubator chamber temperature. This model offers 6" chart of temperature data.



### COA-2013-F 8" Chart Recorder, Temp/Temp, 115/230 VAC, 50/60 Hz

The chart recorder provides an easy-to-read graph of data vs time. It is a reliable, accurate, and stable instrument for on-the-spot written documentation of incubator chamber temperature. This model offers 8" chart of temperature data and comes with 2 remote probes for dual temperature monitoring.




### COA-2014-F 6" Chart Recorder, Temp/RH, 115/230VAC 50/60 Hz

The chart recorder provides an easy-to-read graph of data vs time. It is a reliable, accurate, and stable instrument for on-the-spot written documentation of incubator chamber temperature. This model offers 6" chart of temperature and humidity data.



### 5250001 Voyager Software Kit

Esco Voyager is a PC-based software package developed for the remote monitoring, datalogging and programming / device configuration of Esco controlled environment laboratory equipment. Compatible equipment includes laboratory ovens and incubators, low temperature incubators, CO<sub>2</sub> incubators, and ultra-low temperature freezers.



### COA-2004 (170L) / COA-2042 (240L) 2-UNITS FLOOR STAND

This floor stand allows two units to be stacked without being physically in contact with each other. For the lower unit, it uses roller base for mobility and for easy pull out of the lower unit without the need to remove the upper unit in case of troubleshooting. Floor stand for upper unit also has casters for easy relocation.

# ORDERING INFORMATION

ACCESSORIES	ITEM CODE	DESCRIPTION
COA-1001	5170470	Humidity Display, Factory-installed
COA-1001-F	5170471	Humidity Display, Field-installed Kit
COA-1002	5170472	CO <sub>2</sub> Backup (Tank Switcher), Factory-installed
COA-1002-F	5170473	CO <sub>2</sub> Backup (Tank Switcher), Field-installed
COA-1004	5170474	Reversed Door Swing, Factory-installed
COA-1005	5170475	Analog Outputs, Factory-installed
COA-1005-F	5170476	Analog Outputs, Field-installed
COA-2018-F	5170419	Roller Base (50 L)
COA-2001-F	5170478	Roller Base (170 L)
COA-2019-F	5170420	Roller Base (240 L)
COA-2020-F	5170421	Floor Stand 200 mm (8.0") with Adjustable Feet (50 L)
COA-2002-F	5170479	Floor Stand 200 mm (8.0") with Adjustable Feet (170 L)
COA-2021-F	5170422	Floor Stand 200 mm (8.0") with Adjustable Feet (240 L)
COA-2022-F	5170423	Floor Stand 700 mm (27.6") with Casters (50 L)
COA-2003-F	5170480	Floor Stand 700 mm (27.6") with Casters (170 L)
COA-2023-F	5170424	Floor Stand 700 mm (27.6") with Casters (240 L)
COA-2005-F	5170481	2-Stage Gas Regulator for CO <sub>2</sub> Choose one of the connectors below: 1080588 - CGA 320 Connector (US standard) 1080589 - BP-B534-#8-NT4 Connector (British standard) 1080590 - G5/8-RH Connector (China standard)
COA-2024-F	5170425	Extra Shelf (50 L, Stainless Steel)
COA-2007-F	5170327	Extra Shelf (170 L, Stainless Steel)
COA-2025-F	5170426	Extra Shelf (240 L, Stainless Steel)
COA-2008-F	5170483	Stacking Kit (one set included with every unit purchased)
COA-2010-F	5170329	Electronic CO <sub>2</sub> Analyzer, For CO <sub>2</sub> / Temp Measurement (with Temperature Probe)
COA-2016-F	5170397	Electronic CO <sub>2</sub> + O <sub>2</sub> Analyzer, For CO <sub>2</sub> / O <sub>2</sub> / Temperature Measurement (with Temperature Probe)
COA-2017-F	5170398	Electronic CO <sub>2</sub> + O <sub>2</sub> + RH Analyzer, For CO <sub>2</sub> / O <sub>2</sub> / RH / Temperature Measurement (with Temperature Probe)
COA-2011-F	2170020	IQ / OQ Documentation
COA-2012-F	2170021	6" Chart Recorder, Temp, 115/230 VAC, 50/60 Hz
COA-2013-F	2170022	8" Chart Recorder, Temp/Temp, 115/230 VAC, 50/60 Hz
COA-2014-F	2170023	6" Chart Recorder, Temp/RH, 115/230 VAC, 50/60 Hz
Voyager	5170489	Voyager Software Kit
COA-2004	5250001	2-units Floor Stand Stacking Kit (For 170L)
COA-2042	5170999	2-units Floor Stand Stacking Kit (For 240L)

## After Sales Services

### Parts Availability

Whenever service is needed and parts are required, minimizing downtime is a critical objective. Statistical usage analysis helps Esco to predict parts life, permitting Esco to manage logistics and stage proper inventories around the world. The combination of predictive maintenance, historical data and geospecific proximity assures our customers that parts and labor are available whenever service is scheduled through the local sales organization.

### Registration, Documentation and Instruction

Quality control at Esco extends from research and development through engineering, manufacturing, shipment, delivery and customer feedback. Esco maintains an aggressive program to encourage warranty card registration by mail, email or online submission so that we know where Esco products are located and how they are being used. Rest assured that all information disclosed from warranty registrations will be kept confidential. All Esco products include unique serial numbers for identification. Documentation for all performance tests is archived and maintained for customer reference.

### Online Technical Information

Site preparation instructions are useful before product arrival and installation. Installation and start-up manuals, operation manuals and quick reference guides are available anytime from the Esco resources online.

### NSF International Accreditations

The National Sanitation Foundation (NSF) International is an independent, non-profit organization that provides standards development, product certification, auditing, education and risk management for public health and the environment.

In line with Esco's commitment in providing world class services worldwide, Esco has a large contingent of NSF accredited certifiers which makes Esco not only an **Excellent Standards Company** but also an **Excellent Service Company**, which exemplifies Esco's collective quest of being an **Eternally Successful Company**.

The NSF mark is your assurance that the product complies with all the standard requirements, tested by one of the most respected independent certification organizations in existence today. NSF conducts periodic unannounced inspections and product testing to verify that the product continues to comply with the standard. It is valued by consumers, manufacturers, retailers and regulatory agencies worldwide.

### References and Links

For more information, you can visit Esco at [www.escolifesciences.com](http://www.escolifesciences.com)

Whatever your needs are,  
**we serve it to you  
on a silver platter.**



# ESCO LIFESCIENCES GROUP

42 LOCATIONS IN 21 COUNTRIES ALL OVER THE WORLD



- Global Offices
- Licensee
- Distributors
- Factories
- R&D Centers
- Regional Distribution Centers

*Follow us on social media, download our apps,  
and scan the QR code for more info.*



@EscoLifesciences



@EscoLifesciences



@EscoLifesci



@Esco



@EscoLifesciences



@EscoLifesciences



Esco Lifesciences



Esco Lifesciences

# ESCO

LIFESCIENCES GROUP

Esco Micro Pte. Ltd. • 21 Changi South Street 1 • Singapore 486 777  
Tel +65 6542 0833 • Fax +65 6542 6920 • mail@escolifesciences.com  
www.escolifesciences.com

Esco Technologies, Inc. • 903 Sheehy Drive, Suite F, Horsham, PA 19044, USA  
Tel: +1 215-441-9661 • Fax 484-698-7757  
eti.admin@escolifesciences.com

**Esco Lifesciences Group Offices:** Bangladesh | China | Denmark | Germany | Hong Kong | India | Indonesia | Italy | Japan | Lithuania | Malaysia | Myanmar | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | UK | USA | Vietnam